



Important Customer Notice



NHS Supply Chain

Ref: ICN 2007/0062

Blue Box Medical – Addition to the Epidural Kits and Associated Products Contract

Blue Box Medical will have the following new product available to order via NHS Supply Chain.

Epidural Mini Set

The new 16g Epidural Mini Set will be an addition to the current range of Blue Box Medical Products available via the NHS Supply Chain route allowing a greater degree of choice to our customers.

Components Included within the 16g Epidural Mini Set:

- 1 x 16g x 80mm Tuohy Needle with Wing Attached
- 1 x 19g x 90cm Catheter
- 1 x Standard Loss of Resistance Device Luer Slip
- 1 x Tuohy Borst Adapter
- 1 x Filter 2µ
- 1 x Label



NPC	MPC	Description	UOI
FYBo86	31891-16	Epidural Mini Set 16g 1 x 16g x 80mm Tuohy Needle with Wing attached 1 x 19g x 90cm Catheter 1 x Standard LOR Device Luer Slip 1 x Tuohy Borst Adapter 1 x Filter 2µ 1 x Label	Pack 10

To order this or any other products please contact your local supplies / procurement department or order via your usual method when placing an order with NHS Supply Chain.

Further information

Further information can be obtained from Blue Box Medical, by e-mail enquiries@blueboxmedical.co.uk or contact Marc Buckingham, National Sales Manager by calling 07780-996-883 or Tel: + 44 (0) 2380 669 000, Fax: + 44 (0) 2380 669 002

Visit www.blueboxmedical.co.uk where you can learn more about the products they are able to supply.