

# RS7 Nerve Stimulator

**Special offer valid  
until 31 August 2011**

**BUY 4, GET 1 FREE**  
(In addition to existing  
quantity discounts)

Three stimulation  
options:

Train of Four  
Double Burst  
Tetany 50Hz

Pulse width 200µs

Current options  
between 10mA & 80mA

Handy pocket size



Alarm on incomplete  
circuit

Robust and easy  
to clean

Low battery alarm

Black fabric  
protective case

Probes  
available on request



**BLUE BOX**  
MEDICAL LTD



Call now for a  
demonstration or  
quote, stating ref:  
**RS7 06/11**

Tel: 023 8066 9000 Fax: 023 8066 9002

[enquiries@blueboxmedical.co.uk](mailto:enquiries@blueboxmedical.co.uk)

[www.blueboxmedical.co.uk](http://www.blueboxmedical.co.uk)

Unit 29, New Forest Enterprise Centre

Chapel Lane, Totton, Southampton SO40 9LA

**BLUE BOX MEDICAL**



RS7 06/11