



# AutoFuser® Disposable Silicone Balloon Infuser

- ✓ Continuous infusion for antibiotic, chemotherapy and pain management therapies
- ✓ Easy fill
- ✓ Hard, protective shell
- ✓ Multiple reservoir sizes
- ✓ Self-priming
- ✓ Latex-free
- ✓ Carrying pouch for strap with ambulatory flexibility
- ✓ Optional PFA bolus



**ACEMEDICAL**



# Product Specification



AUTO FUSER <sup>®</sup> DISPOSABLE SILICONE BALLOON INFUSER					
Reference	Sizes		Bolus	Lock-out Time	Units/Box
	Volume	Flow Rate			
AA1505	150ml	5ml/hr	-	-	10
AA1507	150ml	2.5ml/hr	-	-	10
AA1900	100ml	0.5ml/hr	-	-	10
AA1902	100ml	2ml/hr	-	-	10
AA1905	100ml	5ml/hr	-	-	10
AA2002	275ml	2ml/hr	-	-	10
AA2005	275ml	5ml/hr	-	-	10
AA2008	275ml	8ml/hr	-	-	10
AA2010	275ml	10ml/hr	-	-	10
AA2021	275ml	2ml/hr	1ml	15 mins	10
AA2029	275ml	5ml/hr	2ml	30 mins	10
AA4004	400ml	4ml/hr	-	-	10
AA4005	400ml	5ml/hr	-	-	10

This represents a small sample of our comprehensive range of volumes, flow rates and bolus dosage. Please contact us for further details.

