



# Grip-Lok<sup>®</sup> 3300M for General Line Securement

Now available on FP10 Drug Tariff

- ✓ Hypoallergenic, breathable and latex-free skin contact adhesive
- ✓ Provides superior securement for both vertical and horizontal lifting accidental line pulls
- ✓ Soft and flexible fabric design
- ✓ Easy to apply and open
- ✓ Can easily be placed and adjusted with gloved hands
- ✓ Sterile packed to allow application of the device close to the insertion site

Grip-Lok<sup>®</sup> Medium securement device is strong enough to lock lines/tubes in place yet versatile enough for almost any securement situation



| Reference | NPC    | Description              | Tubing Range  | Device Length | Units/ Box |
|-----------|--------|--------------------------|---------------|---------------|------------|
| 3300M     | ELW704 | Grip-Lok Medium 9 - 24fg | 3.0mm - 8.0mm | 8.9cm         | 100 / 1    |



Grip-Lok is a trademark of TIDI Products, LLC

